

PRACTICAL

RUGGED

SMART

# Block Digesters

LABQUEST  
BY BOROSIL®



SAFE

SIMPLE

EFFICIENT

DIGESTION SYSTEMS


*for Labs*

# Block Digesters

SAFE SIMPLE EFFICIENT

DIGESTION SYSTEMS for Labs

- Fast, safe and complete digestion of samples through aluminium blocks
- Intricately control digestion temperature curves for unsupervised digestion
- Monitor process parameters & pre-set sample cycles
- Recommended for industries & Labs that process large numbers of samples in a day, assures complete consistent digestion in all sample tubes



Model Name	KBD062	KBD082	KBD122
Type	Block Digester		
No. of positions	6	8	12
Tube capacity (ml)	250	250	250
Temperature Range	Ambient to 470°C		
Timer Range	1-300 mins		
Temperature Precision	±1°C		
Ramp and Hold	upto 4		
Presets	upto 99		
Auto scrubber ON/OFF	Yes		
PID control	Yes		
Auto OFF after achieving set time	Yes		
Temperature Calibration	Yes, single and two point		
PC & Printer Connectivity	Yes, upgradable		
Display	Graphical LCD		
Body	Stainless steel with control pannel of glass reinforced polymer		
Current Consumption	6A	6A	8A
Dimensions (W x D x H) mm	(414 X 450 X 665)	(416 X 554 X 595)	
Power Consumption	1300W	1300W	1600W
HSN Code	84198990		
Ordering Code	BLFAKBD0620000000	BLFAKBD0820000000	BLFAKBD1220000000

## Cooling Stations



Cooling stations increase throughput by decreasing the cooling time post digestion

**BOROSIL<sup>®</sup>**  
**BOROSIL LIMITED**

1101, Crescenzo, G-Block, Opp. MCA Club,  
 Bandra Kurla Complex, Bandra (E), Mumbai - 400 051, India  
 T : +91 22 6740 6300 F : +91 22 6740 6514  
 E : borosil@borosil.com W : www.borosil.com

Manufactured by :  
**BOROSIL TECHNOLOGIES LTD.**  
 21-24, Plot No. 8, Indialand Global Industrial Park,  
 Hinjewadi Phase - 1, Pune - 411 057

