

Laboratory Bottles & Caps



Made from 3.3 Low Expansion Borosilicate Glass



- Reagent Bottles
- Caps and Pouring Rings
- Aspirator Bottles
- Solution Bottles
- BOD Bottles
- Relative Density Bottles
- Specific Gravity Bottles
- Weighing Bottles

BOTTLES & CAPS

For every laboratory applications

1500
Reagent
Bottle



1501
Reagent
Bottle
Clear



1519
Reagent
Bottle
Amber



SIZES

60 ml	100 ml	125 ml	
250 ml	500 ml	1000 ml	2000 ml

10 ml	25 ml	50 ml	100 ml 150 ml
250 ml	500 ml	1000 ml	2000 ml
3000 ml	5000 ml	10000 ml	20000 ml

10 ml	25 ml	50 ml	100 ml
250 ml	500 ml	1000 ml	2000 ml
3000 ml	5000 ml	10000 ml	20000 ml

Narrow mouth bottle with ground neck and interchangeable glass stopper.
For storage and decanting of acids, alkalis and solvents
Also available in **Amber** colour
cat code 1509

Narrow mouth bottle with screw neck, cap and pouring ring. Provided with retrace code for traceability
Used for preparation and storage of media, solvents and laboratory reagents

10 ml & 25 ml with GL 25 neck
50 ml with GL 32 neck
100ml onwards GL 45 neck

1506
Square
Bottle



1601
Premium
Bottle
Clear



1619
Premium
Bottle
Amber



	500 ml	1000 ml	2000 ml
5000 ml	10000 ml		

			100 ml
250 ml	500 ml	1000 ml	2000 ml
5000 ml	10000 ml	10000 ml	

			100 ml
250 ml	500 ml	1000 ml	2000 ml
5000 ml	10000 ml	10000 ml	

Square bottle with extra wide GL80 neck for easy pouring and filling of liquids as well as powders
Square shape enables stacking of bottles

Also available in **Amber** colour
cat code 1507

Narrow mouth bottle with GL45 neck provided with high temperature resistant cap and silicon liner
Detailed quality certificate available for regulatory compliances
Ideal for storage of formulations, API's, bulk drugs and biopharma products

Choose the right capacity

Correct selection of bottle capacity improves efficiency, space utilization and handling



Bottle Caps and pouring rings

1503

PP Screw Caps
with matching pouring rings



Made of Polypropylene
Autoclavable up to 140°C

1505

PBT Screw Caps
with silicone liner



Made of Polybutylene Terephthalate
Autoclavable up to at 180°C

1513

new

Multiport Caps
with Anti-Twist Mechanism



Made of Polypropylene
Autoclavable up to 140°C

PTFE plugs / caps are available on request

BOROSIL LIMITED

1101, Crescenzo, G-Block, Opp. MCA Club,
Bandra Kurla Complex, Bandra (E), Mumbai - 400 051. India
T : +91 22 6740 6300 F : +91 22 6740 6514
E : borosil@borosil.com W : www.borosil.com

Write to us for export enquiries at eximbor@borosil.com

Laboratory Glassware :
<https://www.borosillab.com>

Laboratory Equipment :
<https://labquest.co>

Filter Papers :
<https://borosilfilters.com>