

PRACTICAL

RUGGED

SMART

Classical Kjeldahl Distillation

LABQUEST
BY BOROSIL®



Image for representation purpose only

SAFE | SIMPLE | EFFICIENT

DISTILLATION SYSTEMS

for Labs

Classical Kjeldahl Digesters

SAFE SIMPLE EFFICIENT

DISTILLATION SYSTEMS for Labs

- Available with complete set of heaters, glassware and accessories
- Classic laboratory heaters for reliable and lasting intensive heating
- Compact for serial heating in batches



CKD065

Model Name	CKD015	CKD065	CKD038*
No. of positions	1	6	3
Flask Capacity	500ml Kjeldahl flask, compatible with 800ml Kjeldahl flask too		800ml Kjeldahl flask
Element Temperature	450 °C		
Temperature accuracy	NA		
Temperature regulators	1	6	3
Ramps/ Flask	4		
Display	LED Display with time & process status		
Timer	1-999 Minutes		
Power consumption	150W	900W	450W
Electrical regulation	0.65A	3.9A	1.9A
HSN Code	84198990		
Ordering Code	BLFACKD01500000000	BLFACKD06500000000	BLFACKD03800000000

* subject to availability

ACCESSORIES

Description- 1 & 6 position	Description- 3 position
Graham condenser 200 mm	Graham condenser 400 mm
Classical Kjeldahl splash head	Classical Kjeldahl splash head
Kjeldahl Flask 500ml	Kjeldahl Flask 800ml

BOROSIL®
Scientific

1101, Crescenzo, G-Block, Opp. MCA Club,
Bandra Kurla Complex, Bandra (E), Mumbai - 400 051, India
T : +91 22 6740 6300 F : +91 22 6740 6514
E : borosil@borosil.com W : www.borosil.com

Manufactured by :
BOROSIL TECHNOLOGIES LTD.
21-24, Plot No. 8, Indialand Global Industrial Park,
Hinjewadi Phase - 1, Pune - 411 057