

PRACTICAL

RUGGED

SMART

Infrared Digesters

LABQUEST
BY BOROSIL®



SAFE

SIMPLE

EFFICIENT

DIGESTION SYSTEMS


for Labs

Infrared Digesters

SAFE SIMPLE EFFICIENT

DIGESTION SYSTEMS for Labs

- Rapid & safe digestion of samples through infrared heaters
- Intricately control digestion temperature curves for unsupervised digestion
- Monitor process parameters & pre-set sample cycles Justified by rapid heating, these digesters find their way in industries that analyse fewer samples.
- Cooling stations increase throughput by decreasing the cooling time post digestion



Model Name	KIR060	KIR020
Type	Infrared Digester	
No. of positions	6	2
Tube capacity	250 ml	
Temperature Controller	Yes	
Element Temperature	700°C	
Time Setting LCD Display	Yes	
Heating Chamber Unit	Stainless Steel	
Current Consumption	12A	4A
Dimensions (W x D x H)	(580 x 390 x 665)mm	(330 X 360 X 665)mm
Power Consumption	2400W	800W
HSN Code	84190000	
Ordering Code	BLFAKIR06000000000	BLFAKIR02000000000
Prices / Piece	On request	

Cooling Stations



THW020



THW062

Cooling stations increase throughput by decreasing the cooling time post digestion

BOROSIL
BOROSIL LIMITED

1101, Crescenzo, G-Block, Opp. MCA Club,
Bandra Kurla Complex, Bandra (E), Mumbai - 400 051, India
T : +91 22 6740 6300 F : +91 22 6740 6514
E : borosil@borosil.com W : www.borosil.com

Manufactured by :

BOROSIL TECHNOLOGIES LTD.

21-24, Plot No. 8, Indialand Global Industrial Park,
Hinjewadi Phase - 1, Pune - 411 057