

Multi position stirrers with / without heating

5 P / 10 P / 15 P

Multi position stirrers with slip-less speed regulation, homogenous temperature distribution and brushless motor for maintenance free operation.

The closed design allows for easy cleaning, safety and protection against spillage.





MHPS 5 P



MHPS 10 P



MHPS 15 P

Technical Information

Specification	MHPS 5P	MHPS 10P	MHPS 15P
Stirring positions	5	10	15
Speed (rpm)	1500	1500	1500
Display	digital	digital	digital
Incremental speed (from 15 rpm)	10	10	10
Max capacity in lit	2	4	6
Temp range (from ambient)	120°C	120°C	120°C
Dimension (W X D X H) in mm	562 x 136 x 75	565 x 196 x 75	568 x 282 x 76
Set up plate in mm	456 x 118	460 x 178	460 x 264
Pulse mode	Yes	Yes	Yes
Product code (with heating)	100MHPS0515000	100MHPS1015000	100MHPS1515000
Product code (without heating)	100MS000515000	100MS001015000	100MS001515000

BOROSIL®

Pioneers in speciality glass

www.borosil.com

Borosil Glass Works Ltd.

Corporate Office

1101, Crescenzo, G-Block, Opp. MCA Club, Bandra Kurla Complex, Bandra (E), Mumbai - 400051, India
Tel : +91 22 6740 6300 • Fax : +91 22 6740 6514 • Email : borosil@borosil.com

For export enquiries write to us eximbor@borosil.com