



Still stuck with
the old method

for transferring HPLC solvents ?

BOROSIL[®]

PROUDLY MADE IN
india

is right here with a solution

Multi-Port HPLC Cap

with Anti-twist mechanism



*A smart solution
for transfer of
HPLC mobile
phase solvents*

How **BOROSIL** Multi-Port Caps keep your lab Environment safe ?



Ensures a safe and closed loop transfer of HPLC solvents

Prevents spillage and leakage of your expensive solvents during transfer

Prevents release of harmful solvent vapours in the laboratory

1513 - Multiport HPLC Cap

Unique air valve cum filter

0.1µm filter prevents contamination of solvents with dust particles and the release of solvent vapours thereby maintaining the solvent concentration

PP fittings for different tubing ODs

Supplied with three connectors – for 3.2 mm and 1.6 mm OD tubing alongwith one end plug

Polypropylene Cap

Suitable for all types of GL 45 bottles and autoclavable at 121°C

Anti-twist mechanism

Freely rotatable PTFE manifold prevents twisting of the tubes while opening and closing the cap



Product Code	No of Ports	Neck	Quantity per Pack
1513C04	4	GL45	4 Nos
1513C02	2	GL45	4 Nos
1513C01	1	GL45	4 Nos

Ideal for use with

Reagent Bottle

and HPLC Solvent Bottle



1501 Reagent Bottle with GL 45 Neck

Product Code	Capacity ml	Approx O.D. x Height*	Thread Specification	Quantity Per Case
1501016	100	56 x 100	GL 45	10
1501018	150**#	60 x 120	GL 45	10
1501021	250	70 x 138	GL 45	10
1501024	500	86 x 177	GL 45	10
1501029	1000	101 x 225	GL 45	10
1501030	2000	136 x 260	GL 45	10
1501031	3000**	158 x 295	GL 45	2
1501033	5000***	181 x 320	GL 45	1
1501038	10000***	235 x 405	GL 45	1
1501040	20000***	296 x 505	GL 45	1

- *Height indicated is of Bottle only
- #Without Pouring Ring
- **Not covered in any standard
- ***As per Borosil Dimension

Multi-Port HPLC Cap

with Anti-twist mechanism



*A smart solution
for transfer of
HPLC mobile
phase solvents*

BOROSIL®

enabling the future of science

BOROSIL LIMITED

1101, Crescenzo, G-Block, Opp. MCA Club,
Bandra Kurla Complex, Bandra (E), Mumbai - 400 051. India
T : +91 22 6740 6300 F : +91 22 6740 6514
E : borosil@borosil.com W : www.borosil.com

Write to us for export enquiries at eximbor@borosil.com