

Still using
non-certified vials for
storage of expensive
working standards ?

BOROSIL[®]

PROUDLY MADE IN
india

is right here with a solution

Certified Sampling and Storage Vials

*An ideal container for storage of
working standard solutions*



Made from 5.1 Low Expansion
Borosilicate Glass



How Borosil Storage Vials different from others ?



Provided with a Quality Certificate

Packed in ISO "Class 8" clean room facility

Complies glass composition and Spectral Transmission test as per USP chapter 660/ USP 40 and European pharmacopoeia 9.0

9913 - Sampling & Storage Vials



PP Screw Cap with PTFE liner for leak proof applications

Amber colour to protect light sensitive chemicals

Write-on Patch for easy labeling and identification

Made from USP Type I, 5.1 expansion borosilicate glass (Glass Type as per ASTM E-438)

Product Code	Capacity ml	Approx Mouth ID mm	Dimensions O.D. x Length mm	Case lot (nos)
9913003	3	10	16 x 36	1050
9913005	5	10	16 x 55	1050
9913006	10	16.5	24 x 50	600
9913007	15	16.5	24 x 58	600
9913008	20	16.5	24 x 70	600
9913010	30	16.5	24 x 95	600

*Also available in Clear Type



BOROSIL LIMITED

1101, Crescenzo, G-Block, Opp. MCA Club, Bandra Kurla Complex, Bandra (E), Mumbai - 400 051, India
 T : +91 22 6740 6300 F : +91 22 6740 6514
 E : borosil@borosil.com W : www.borosil.com

MUMBAI

T : 022 6740 6400
 F : 022 6740 6444
 E : mumbai@borosil.com

DELHI

T : 011 41505893
 E : delhi@borosil.com

CHENNAI

T : 044 2822 6012 / 13
 F : 044 2822 6014
 E : chennai@borosil.com

KOLKATA

T : 033 2229 9166
 F : 033 2226 2045
 E : calcutta@borosil.com

Write to us for export enquiries at eximbor@borosil.com

Laboratory Glassware : <https://www.borosillab.com>