


BE VOCAL ABOUT LOCAL

# Soxhlet Apparatus



PRACTICAL

RUGGED

SMART

# The Only Soxhlet Apparatus With ACTUAL TIME SETTING AND ENERGY PERCENTAGE DISPLAY

Carry out soxhlet extraction up to 6 samples in parallel with the series Gen 2 heating mantles. The functional, robust construction keeps your running costs low and at the same time minimizes the risk of failure. The compact design ensures work safety and ergonomic working

## SAFE

Inbuilt safety control



## INTELLIGENT

LED display of time, safe temperature and energy percentage




## PRECISE

Inbuilt Precise temperature control upto +/- 1 °C



## START YOUR SOXHLET AND WALK AWAY!


# The smart SOX series of soxhlet extractors

- 
- Extractor knob for solvent recovery
  - Digital timer with 1 to 999 min time setting.
  - Digital energy regulator with 10 step power regulation.
  - Safe temperature setting to prevent accidents.
  - Smart stand for condenser placement facilitating ease of cleaning
  - Auto heater cut off with visual & audio error indicator to user
  - Chemical and heat resistant body enhancing life of unit
  - Built in calibration for error setting.
  - Large character digital display for better visibility.
  - Rotary knob for precise setting of time, safe temperature and energy percentage regulation

# Technical Specifications



SOX012  
1 position soxhlet



SOX062  
6 position soxhlet

Model Name	SOX012	SOX015	SOX042 *	SOX045*	SOX062	SOX065
Position	1	1	4	4	6	6
Extractor capacity (ml)	100	200	100	200	100	200
Flask Capacity (ml)	250	500	250	500	250	500
Temperature Accuracy	±1°C					
Mantle Temperature	Upto to 320°C					
Safe temperature setting	Yes					
Knob for Extractor	Yes					
Display Timer	1-999 min & Continuous					
Dimension (W x D x H)mm	200 x 300 x 800		750 x 300 x 800		1050 x 300 x 800	
Weight (Kgs)	7		12		15	
Power (W)	400		1600		2400	
HSN Code	84194090					
Ordering code	BLFRSOX0120 0000000	BLFRSOX0150 0000000	BLFRSOX0420 0000000	BLFRSOX0450 0000000	BLFRSOX0620 0000000	BLFRSOX0650 0000000

\*Till stock lasts

**BOROSIL®**  
**BOROSIL LIMITED**

1101, Crescenzo, G-Block, Opp. MCA Club,  
Bandra Kurla Complex, Bandra (E), Mumbai - 400 051, India  
T : +91 22 6740 6300 F : +91 22 6740 6514  
E : borosil@borosil.com W : www.borosil.com

Manufactured by :

**BOROSIL TECHNOLOGIES LTD.**  
21-24, Plot No. 8, Indialand Global Industrial Park,  
Hinjewadi Phase - 1, Pune - 411 057

LABQUEST  
BOROSIL

BOROSIL  
TECHNOLOGIES

