

GLASSWARE DRYER





Easy to Use



Adjustable Cycle Timer



Safety



Low Maintenance

Labquest Glass Dryers can be used as benchtop models or wall hanging in your laboratory washing area. It offers drying of laboratory articles in 10 to 15 minutes. Useful for Quick Drying Viz. GCMS Glasswares, Kjeldahl Digestion Tubes, Routine work involving Assorted Glassware, General Purpose Drying.

Key Features:

- Light Weight
- Inter-Changeable PEGS
- 0-60 mins Cycle Operation
- Versatile Glassware Compatibility
- Replaceable Air Filter

Technical Specs

Feature/ Model	GDS020
Number of PEGS	21
PEG Sizes	
16 mm	8 Nos.
10 mm	8 Nos.
7mm	5 Nos.
Temperature Range	Ambient to 70°C
Time Control	Yes
Heater Selector	Manual
Drying Power	650 W
Dimensions in mm (L x B x H)	440 x 170 x 550
Voltage	240 V
Ordering Information	BLFGGDS02000000000
HSN Code	84798200

