

BE VOCAL ABOUT LOCAL

# Magnetic Hot Plate Stirrer

**550 PREMIER**



The nano crystalline ceramic glass-top provides excellent resistance to thermal shock and corrosive liquids. The elevated design, protects the unit from spillages. Best in class capacity and speed of 20 lit and 2200 rpm. The unit comes with slip detection that can recover the stir bar in the event of it being lost outside the magnetic field. The target temperature can be achieved in rapid or gradual or accurate modes depending upon the application. The unit can also be programmed to meet custom temperature and time conditions.

## Magnetic Hot Plate Stirrer

**550 PREMIER**



### Specifications :

Ordering Information	100MHPS010550P
Max Capacity	20L
Max Speed	200-2200 rpm
Set-up Plate	Glass Ceramic
Setup Plate Dimensions	180 x 180
Max Plate Temp	550
Max Liquid Temp	280
Stirring Positions	1
Timer	1-999min & Infinite
Digital Display	Set & Actual Temp / RPM / Time
Safe Temp	Yes
Hot Indication	Above 50°
Protection Class	IP 21
Dimension in mm	330 x 220 x 115
Stir Bar	25 mm
Pulse Mode	Yes
External Probe	PT 1000
HSN Code	84198990

**BOROSIL**  
**BOROSIL LIMITED**

1101, Crescenzo, G-Block, Opp. MCA Club,  
Bandra Kurla Complex, Bandra (E), Mumbai - 400 051. India  
T : +91 22 6740 6300 F : +91 22 6740 6514  
E : borosil@borosil.com W : www.borosil.com

