

ON DEMAND
Quartz Water Distillation Units



PRACTICAL

RUGGED

SMART

ALL GLASS

Quartz Double Distillation Unit

ON DEMAND WATER DISTILLATION UNITS

for labs



Save Water



Power & Water



Safety



Low Maintenance

- On Demand Distillation - get the exact volume of distilled water without intervention or monitoring.
- Cleaning guide for efficient operation.
- Lower power and water consumption.
- Digital output display.

Purified Water Specs - Typical

Feature/ Model	QUARTZ DDU
Conductivity (S/ cm)	< 1 x 10 ⁻⁶
Biological Activity	Pyrogen Free

Technical Specs

Feature/ Model	QUARTZ DDU
Dist. water output cap (Approx) (Ltr/hr)	2.5 LPH
Minimum cooling water requirement (Ltr/min)	1.5
Power consumption (kw)	5.5
Distillation Temp (°C)	65 - 75°C
Voltage (V)	230 - 250V
Heater	Quartz
Boiler	Quartz
Condenser	Quartz
Dimensions in mm (W x D x H)	650 x 630 x 520
HSN Code	84194090
Ordering Information	BLFRQDD25000000000

* Inlet water TDS < 100 ppm

Note: Because of difficulties associated with measurement of the pH value of high purity water and the doubtful significance values obtained, limits for the pH value of grade 1 & grade 2 water have not been specified.

BOROSIL®
BOROSIL LIMITED

1101, Crescenzo, G-Block, Opp. MCA Club,
Bandra Kurla Complex, Bandra (E), Mumbai - 400 051, India
T : +91 22 6740 6300 F : +91 22 6740 6514
E : borosil@borosil.com W : www.borosil.com

Manufactured By :

BOROSIL TECHNOLOGIES LTD.

Plot No. 7, Sr. No. 234, 235, & 245,
Indialand Global Industrial Park, Hinjewadi, Phase 1,
Pune - 411057, Maharashtra, India.

LABQUEST
BOROSIL

BOROSIL
TECHNOLOGIES

