

enabling the future of science

BOROSIL[®]
Scientific

Borosil Scientific Limited



PROTECT
your HPLC system
with demanding

HPLC
Inlet Filters

Borosil HPLC Inlet filters are designed to provide optimal performance and prevent contamination. These filters act as a shield, protecting the system and prolonging its lifespan.

PROTECT

your HPLC system
with demanding

HPLC Inlet Filters

- Made of high-quality 316 L sintered stainless steel
- Controlled porosity to selectively retain particles
- Excellent chemicals and corrosion resistance
- Compatible with most HPLC systems & tubing sizes
- Bubble point tested to verify pore size and integrity for every lot of production



Universal Slip on 1/8



Push Fit PTFE Solvent Filter



Step Tube Fit Solvent Filter



Note: Use metal filters only with compatible solvents and avoid contact with corrosive chemicals.

Product Code	Description & Dimensions	Filter rating	Case Lot Units / Pack
SIFBSSW1001	SS316L-SS 316-SLIP ON-1/8" solvent filter, single step, filter rating: 10 micron, size: 126425, MOC: SS 316L / SS 316L, 1/8" Slip-on	10 um	100
SIFBSSS1002	SS316L-SS316 push fit PTFE, solvent filter, single step, filter rating 10 micron, size 126-425, MOC-SS316L, Push Fit PTFE	10 um	100
SIFBSSSU1003	SS316L-SS316 Step Tube Fit, solvent filter, single step, filter rating 10 micron, size 126-425, MOC-SS316L, Step Tube Fit	10 um	100

* 2micron Filter rating also available on request , MoQ exists

All filters are for intended application use only

BOROSIL®
Scientific
Borosil Scientific Limited

CIN : U74999MH1991PLC061851

Registered & Corporate Office :
1101, Crescenzo, G-Block, Opp. MCA Club, Bandra Kurla Complex,
Bandra (E), Mumbai - 400 051, India.
T +91 22 6740 6300 • F +91 22 6740 6514
E bsl@borosil.com • W www.borosilscientific.com
Write to us for export enquiries at bsl.exports@borosil.com

MUMBAI

mumbaisales@borosil.com
+91 22 6740 6400
+91 22 6740 6444

DELHI

delhi@borosil.com
+91 11 4150 5893

CHENNAI

chennai@borosil.com
+91 44 2822 6012 / 13
+91 44 2822 6014

KOLKATA

calcutta@borosil.com
+91 33 2229 9166