


Rapid Fat Analyser



Rapid fat extractor with benchmarking features of auto solvent Recovery, Boiling & Rinsing

Carry out fat extraction up to 6 samples in parallel with the series. The functional, robust construction keeps your running costs low and at the same time minimizes the risk of failure. The compact design ensures work safety and ergonomic working

SAFE

Inbuilt safety control



INTELLIGENT

Touch screen display to monitor and control all functions




PRECISE

User defined setting of rinsing, boiling and recovery time



START YOUR UNIT AND WALK AWAY!

The smart rapid fat extractors



Technical details

Model Name	FAS060
No of positions	6 No's
Sample size	0.5 gm to 10 gm
Recovery of Solvent	80-90%
Limit of detection	0.1% to 100 %
Heating profile	Energy regulator 1 to 10
Display	7" touch screen
Water flow sensor	Yes
Beaker solvent capacity	150 ml
Condenser	Seal Viton (FKM)
Chiller water temp	10-15 degree
Heater temp range	280 degree Max
Recovery bottle capacity	2 liter
Flow cut off sensor	yes
Vacuum pump flow rate	380 L/hr
Power supply	230 V 50 – 60 Hz
Current Consumption	5.2A
Power consumption	1200 W
Dimensions (LxBxH) mm	840 x 430 x 600
HS Code	84198990
Ordering code	BLFRFAS060000000000
Price/ Piece	On request