

Ultrasonic Baths

Ultrasonic Bath offers gentle effective cleaning action ensuring that no contaminants get carryover from the previous process.

Sonicator Bath

Ultrasonic Baths - Use sound waves and liquid to clean small objects. Borosil Labquest digital series of Ultrasonic cleaners come in four different sizes. These Stainless Steel cleaners feature a bright LED display for time and temperature and easy push-button control along with degs functions.



Sonicator Bath - Features

- Temperature control from Ambient - upto 80 °C

Small Footprint, Supplied with Lid
- Display of target temp, actual temp & time in a Microprocessor Controlled 7 segment Digital Display

Time setting upto 99 mins or in continous mode
- Degas Mode - better Cleaning

Seamless tanks ensuring better cleaning

Drain Tap for 5L & above models



Technical Data & Ordering Information

Model Name	UCB025	UCB050	UCB100	UCB200
Capacity	2L	5L	10L	20L
Working Temperature Range	Ambient +5 °C to 80 °C			
Frequency	40kHz			
Time Setting	1- 99 mins	1- 99 mins	1- 99 mins	1- 99 mins
Interior Dimension (L, B, H) in mm	150 x 140 x 100	300 x 155 x 150	300 x 240 x 150	500 x 300 x 150
Exterior Dimension (L, B, H) in mm	190 x 170 x 220	330 x 180 x 330	330 x 270 x 330	550 x 330 x 330
Ultrasonic Power	50 W	120W	240W	500 W
Heater Power	100W	100W	500 W	500 W
Ordering Code	UCB02500000	UCB05000000	UCB10000000	UCB20000000

HSN Code : 84798999